## **HyLED 9500 Mobile**

## LED Surgical Lights

The materials and information provided via this brochure are distributed internationally, but not all the products or services found on this brochure are necessarily available in your country or area. Please contact local office or agents to make sure of the availability of the relevant products or services.

Naddress: Rua Joaquim Floriano, 488, floor Intermediário - cj 02, Itaim Bibi, São Paulo - SP, CEP 04534-002, Brazil Tel: (55-11) 3078-8026

Address: 4250 Kingsway, Suite 206, Burnaby, BC, Canada V5H 4T7 Tel: (1-604) 451-1199

**Egypt** Address: 24th floor,34 EL-Obour Buildings, Salah Salem Street, Helioplis, Cairo, Egypt Tel: (20-2) 2260-9728

France Address: Europarc Créteil Cedex 123, Chemin des Bassins 94000

Tel: (33-1) 4513-9145

Address: Zwischen den Bächen 4, 64625 Bensheim Tel: (49-6251) 17524-0

Address: Unit No. 401/402, NDM-1, Netaji Subhash PlaceWazirpur District Centre, New Delhi 110 034 India Tel: (91-11) 4923-0000

Indonesia Address: Menara Dea 1, Suite 301, Kav E4 No.3/1, Mega Kuningan, Jakarta Selatan 12950, Indonesia

### Address: Via Girardo Patecchio, 4 - 20141 Milano, Italy Tel: (39-02) 57402-4444

Address: Insurgentes Sur # 1787 Piso 2, Col. Guadalupe Inn, Mexico D.F. 01020

Tel: (52-55) 5661-9450

Address: Saturnusstraat 17, 2132HB Hoofddorp, The Netherlands Tel: (31-23) 574-9656

Address: Kondratyuka Str. 3 129515, Moscow Tel: (7-499) 553-6036

### Address : Avenida Manoteras 38, Bloque B, Bajo B008. Madrid

28050, Spain Tel: (34-91)392-3754

Address: General Ali Rıza Gürcan Cad. Metropol Center No:31 A Blok K.13 D.52 Merter/Istanbul, Turkey

Address: 3 Percy Road, St. John's Park, Huntingdon, Cambs, Pe29 6SZ, United Kingdom Tel: (44-1480) 416-840

DISTRIBUTOR:

Address: 800 MacArthur Blvd., Mahwah, NJ 07430-0619, USA



# mindray

Mindray is listed on the NYSE under the symbol "MR" Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China Tel: +86 755 26582888 Fax: +86 755 26582680 E-mail: intl-market@mindray.com Website: www.mindray.com





# **HyLED 9500 Mobile**

## LED Surgical Lights

HyLED 9500 Mobile offers long service time, excellent LED illumination, superior flexibility and manoeuvrability, which caters for emergency, minor surgery, examination and treatment rooms, and operating rooms.





High quality LED bulb creates a homogenous light beam



Really thin, the lateral profile Of the light head



Cross shaped design for better shadow dilution and laminar flow compatibility



Simplicity design of the handle provides superior maneuverability



The ultra-agile mobile base provides maximum flexibility, with reliable battery power supply up to 1.5 hours



Four castors all with locking system provide stable base for the light

### **Technical Specifications**

	HyLED 9500M
Central illuminance (at 1m distance)	120,000 lux
Light field diameter(at 1m distance)	195-330mm
Depth of illumination (L1+L2)	1,200mm
Color temperature	4,350K
Color rendering index	95±3
Shadow dilution with tube	100%
Shadow dilution with one mask	65%
Shadow dilution with tube and one mask	60%
Shadow dilution with two masks	55%
Shadow dilution with tube and two masks	50%
Radiant energy	≤4.5 mW/m²*lux
Ambient illumination	9,800 lux
Temperature rise-area around surgeon's head	≤1℃
Temperature rise-operative field	≤1℃
Service life of bulb	≥40,000 hours
Power consumption	150VA
Power supply	100VAC~240VAC,50/60Hz
Battery supply duration	1.5 hours(full power) 3 hours(normal state)

